

## **ID 162VB**



- Pre-sterilized, disposable single piece implanter
- Biocompatible glass encapsulation
- Non-contact animal ID
- ISO 11784/85 FDX-B compliant
- 6x adhesive label featuring ID number in numeric and barcode 39 format
- Implanter provides tactile "click" ensuring correct positioning of the tag in the animal.

## **ID 162VB**

**Animal Transponder** in VetPlant<sup>TM</sup> Implanter FDX-B ISO 11784/85

The ID-162VB Animal Transponder in VetPlant Implanter is a miniature transponder, encapsulated in biocompatible glass for use in animal applications. Each ID-162VB transponder is ready for application, in its own disposable needle and 6 ea. adhesive labels featuring the transponder's ID number in barcoded format and alpha-numeric characters.

The ISONORM system features FDX-B Transponders compliant with ISO 11784 and ISO 11785 for use in animal identification applications which do not require ID code security.

The VetPlant<sup>TM</sup> Implanter provides tactile feedback when the transponder is about to move into the eye of the needle, which helps preclude premature extrusion. The Implanter provides a second tactile "click" when the transponder is fully extruded, so the operator knows "it's ok to pull out the needle." This feature precludes the transponder from coming back out with the needle. Once fully extruded, the push pin also locks in place at the point of maximum extrusion so it doesn't suck the transponder back into the needle.

## TRANSPONDER

DIMENSIONS

11.5 mm (0.44 in.) Length: Diameter: 2.12 mm (0.083 in.)

WEIGHT

 $0.095g \pm 0.02g$ 

**IDENTIFICATION CODE** 

128 bits factory pre-programmed

**PROTOCOL** 

ISO 11784/ISO 11785 FDX-B

**ENCAPSULATION** 

Biocompatible glass

SCAN ANGLE

Spherical

STORAGE TEMPERATURE (tested exposure 1 hour)

 $-40^{\circ}$ C to  $+90^{\circ}$ C ( $-40^{\circ}$ F to  $194^{\circ}$ F)

OPERATING TEMPERATURE

 $-25^{\circ}$ C to  $+85^{\circ}$ C ( $-13^{\circ}$ F to  $185^{\circ}$ F)

OPERATING FREQUENCY

134.2 kHz

**ENVIRONMENT** 

**IP68** 

MECHANICAL PRESSURE

Axial max. 10N

Radial max. 40N

**VIBRATION** 

Random -20 to 2,000 Hz/10 g./1 hr. p. axis Sine -20 to 2,000 Hz/5 g/1 hr. p. axis

NEEDLE

**DIMENSIONS** 

Length:

32 mm (1.25 in.) 2.6 mm (0.10 in.)

Outside dia.:

SYRINGE ASSEMBLY INC. NEEDLE

138.5 mm (5.45 in.) Length Max width: 45 mm (1.77 in.)